Ursprungssprache

DE

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
## EN

### Safety instructions

- Read the operating instructions in full before starting up and follow the safety instructions.**

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Refer to the operating instructions for the drive.
- Abrasion of the milling unit or the rotating accessories can get into the medium you are working on.
- Do not use any damaged tubes e.g. with hairline cracks or splits.
- Excessive sizes and excessive hardness can damage the **MMT 40/100** (see **"Technical data"**).
- Twist the milling tube on securely until it clicks home audibly.
- Take care not to fill the vessel above the mark.
- Ensure that the tube is firmly attached to the bayonet lock connector on the drive unit prior to operating the device.
- After the milling operation, the **MMT 40/100** may be hot. Please either allow it to cool down or wear appropriate personal protective equipment before detaching it from the drive.
- The processing time is limited to three minutes. The **MMT 40/100** must be allowed to cool down for a few minutes before processing the sample further in the same mill.
- The temperature of the material must not exceed 125 °C. ***Caution!*** A sample of dry maize may reach this temperature just by milling for three minutes at maximum speed.
- The use of liquid nitrogen for directly cooling the sample in the sample vessel is not permitted.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Do **not** use the tubes in explosive atmospheres, with hazardous substances or under water.
- Do not process materials that are radioactive, flammable or explosive, or materials that react chemically with each other releasing large amounts of energy.

## Product information

- MMT 40/100** milling tubes are for use in conjunction with the drive **Tube Mill control/Tube Mill 100 control**, for milling dry samples up to 5 Mohs hardness. The sample pieces must not be larger than 1 cm³. If the sample pieces are larger than this they must be pre-milled before placing them in the **MMT 40/100**. The harder the sample pieces, the smaller they must be. If for instance the hardness is 5 Mohs the edge length must be no longer than 0.5 cm (corresponding to a piece size of approx. 0.125 cm³). The tubes can be used for preparing samples for e.g. quality control in food industry or pharmacy. Commodities suitable for processing include grains, sugar, coffee, tea and tablet etc..
- Dry ice can be used to extend the application range to include moist and resilient samples such as fruits, leaves, jelly babies etc..

- Delivery scope:**

- Multi-use milling tube <b>MMT 40/100</b>	- 5 pieces beaters	- 5 pieces couplings
- 25 pieces sealing rings	- Torque handle	- Tube fixture
- Operating instructions.	Accessories see <b>www.ika.com</b>	

## Application instruction

- Store the tubes in a dry place, in their packaging and at room temperature.
- Protect the tubes against ultraviolet light.
- Do not fill the **MMT 40/100** above the maximum mark on the sample vessel. Milling less material may give better results (eventual finer particle size).
- If the samples are resilient or contain residual moisture, add a teaspoonful of dry ice (dry ice snow, not compacted dry ice) directly into the sample vessel.
- To close the **MMT 40/100**, twist the coupling cover clockwise until it clicks home audibly, see **Fig. 1**.
- The tube is attached by twisting it clockwise on to the bayonet coupling of the drive. The tube is correctly attached when it has clicked home audibly, see **Fig. 2**.
- When the coupling cover has been correctly screwed on to the tube and the tube is correctly attached to the drive, one of the grip cams on the sample vessel will point directly to the rear, see **Fig. 3**.
- Longer processing times will cause high temperatures in the milled material, for instance a 40 ml sample of maize milled for three minutes at maximum speed can reach a temperature of approx. 125 °C.
- If the sample is temperature-sensitive, run the mill for only a short time. Cooling the sample with dry ice snow can also prevent it overheating.
- Likewise the optimal process duration and rotating frequency must be determined by attempts.
- Before you open the **MMT 40/100** it is advisable to tap lightly on the table so that all the milled material falls into the sample vessel.
- To open the **MMT 40/100**, release cam **(A)** and twist the coupling cover anticlockwise **(B)**. Open the tube always up side down (see **Fig.4**). ***Attention:*** Don't open the tube in the operation position.

## Replacing the beater, the coupling and sealing ring

- When the beater or coupling is used for long time , remove them as shown in **Fig. 5** for cleaning and replacing the sealing ring:
  - ⇨ Place the coupling cover into the tube fixture.
  - ⇨ Align the notch in the torque handle with the beater and turn clockwise (left-hand thread) to unlock.
  - ⇨ The beater can now be removed from the coupling cover, leaving the sealing ring and coupling also ready for removal.

- Assemble the new coupling, sealing ring and beater as show in **Fig. 6**:
  - ⇨ Place the coupling and coupling cover into the tube fixture.
  - ⇨ Lay the sealing ring in the recess in the coupling cover and press down into position.
  - ⇨ Place the beater in the coupling recess and lock in place by turning the torque wrench anti-clockwise (left-hand thread).
  - ⇨ The beater is locked in position when a **"click"** can be heard.

***Note:*** The sealing ring must be replaced after each grinding process. The sealing ring is for single use. A cleaning process isn't possible for the sealing ring.

Technical data		
	<b>MMT 40</b>	<b>MMT 100</b>
Sample vessel and coupling cover	Polyamide (PA)	
Inner vessel	Stainless steel (AISI 316L)	
Beater	Stainless steel (AISI 301)	
Sealing ring	Polyethylene (PE)	
Coupling	Polyoxymethylene (POM)	
Max. filling volume	40 ml	100 ml
Max. circumferential speed	65 m/s	
Max. material size	1cm edge length	
Max. material hardness	5 Mohs	
Max. processing time	3 minutes	
Max. temperature	125 <span> </span> °C	
Suitable for	Tube Mill control /100 control	Tube Mill 100 control
Suitable for dishwasher cleaning	all parts except for the sealing ring	
FDA conformity	Yes	
Sterilization (121 <span> </span> °C, 2 bar)	all parts except for the sealing ring	

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