Procedures for	Gel	Preparation
with SequaGe	I XR®	

SequaGel XR (EC-842) is a convenient two bottle gel system. SequaGel XR Monomer Solution contains urea as well as acrylamide and acrylamide derivatives (in a proprietary ratio) dissolved in deionized, distilled water. SequaGel Complete Buffer Solution contains 5X TBE and TEMED in deionized, distilled water. Store solutions tightly capped in a dark area at room temperature (20°C). If SequaGel Monomer Solution is refrigerated, urea may precipitate. This urea will redissolve when the solution returns to room temperature. Acrylamide has been found to be neurotoxic. Protective gloves should be worn while handling these products. If accidental exposure occurs, contact a physician immediately.

Mix SequaGel XR Monomer Solution and Buffer Add appropriate volumes of SequaGel XR Monomer Solution and UreaGel Complete Buffer to a thick-walled Erlenmeyer flask (see Table 1). If desired, the solution may be degassed by stirring under vacuum for two minutes. Bring to room temperature before polymerization.

Table 1: Volumes of SequaGel XR Monomer Solution and SequaGel Complete Buffer to prepare 100mL gel solution			
SequaGel Monomer Solution	SequaGel Complete Buffer	10% Ammonium Persulfate	
80 mL	20 mL	800 µL	

Add APS and Cast Gel

Suggestions for Best Results

Add FRESHLY PREPARED 10% ammonium persulfate (see Table 1). Swirl gently to mix, and cast the gel. Insert the comb and allow to polymerize one to two hours. NOTE: After two hours of polymerization wrap each end of the

 Clean glass plates thoroughly. Rinse with ethanol and wipe dry. Apply Glass Free (Cat.# EC-621) to one plate to ensure release after electrophoresis.

- Prerun the gel for 15-30 minutes before loading the samples. The gel temperature should be between 45-50°C.
- After the completion of the run, allow the plates to cool 10-15 minutes before separation.

Order No. EC-860

gel cassette with clear plastic wrap. This is important to keep the ends of the gel from drying and to maintain sample well integrity. Appropriately wrapped gels may be stored for up to 48 hours.

Table 2: Tracking Dye Migration in SequaGel XR			
Bromophenol (nucleotide	Blue s)	Xylene Cyanole (nucleotides)	
50		210	
			_

Use the table above to monitor electrophoresis progress by means of the co-migration of the bromophenol blue and xylene cyanole tracking dyes with your samples. When doing multiple loads, the next load should be added when the bromophenol blue is 3-4 cm from the bottom of the gel.

Available from National Diagnostics for Nucleic Acid Separation and Analysis

10X TBE

Strict quality control from raw materials to finished product

- Saves time
- 18 MegOhm water
- 0.2 micron microfiltration
- Good price
- Excellent, focused technical support



National Diagnostics' is the purest and most stable on the market. 10X TBE contains 0.89M Tris Borate pH 8.3 and 20mM Na₂EDTA.

3H 14C 3H 14C

Autofluor^{T.} Order No. LS-315

High Resolution Autoradiographic Image Intensifier

- Rapid enhancement of low energy betaemitters such as ³H, ¹⁴C, and ³⁵S
- For polyacrylamide gels, paper chromatography, and TLC plates
- Water based, odorless, contains no DMSO

PPO-DMSO Autofluor

CONTACT NATIONAL DIAGNOSTICS

EUROPE: 0115 982 1111