

R-Luc mRNA

(mRNA encoding Renilla Luciferase protein)


Description

Ready-to-use stabilized R-Luc mRNA
 Concentration: 1.0 mg/mL in 1 mM Sodium Citrate, pH 6.4.
 mRNA length: 1150 nt. MW **MRNA20**= 373000 g/mol; **MRNA21**= 377650 g/mol; **MRNA7**= 375325 g/mol.

R-Luc mRNAs have been designed to produce high expression level of Renilla Luciferase protein. OZB mRNAs are produced by *in vitro* transcription. mRNAs are stabilized at the 5' end by modified nucleotides capping (Cap1) and contain a poly(A) tail at the 3' end. Sequences have been optimized to yield improved stability and performance. R-Luc mRNA #**MRNA20** does not bear any additional nucleotide modifications while #**MRNA21** is modified with 5-methoxyuridine (5moU), #**MRNA7** is modified with N1-methyl-pseudouridine (N1-mψ) to reduce innate immune response.


(ref# **MRNA20**):

Mature mRNA (unmodified nucleotides) with cap1 and polyA tail




(ref# **MRNA21**):

Mature mRNA (fully modified moU) with cap1 and polyA tail



(ref# **MRNA7**):

Mature mRNA (fully modified N1-mψ) with cap1 and polyA tail



Applications

R-Luc mRNAs can be used as control of transfection efficiency as a reporter gene. R-Luc mRNAs resemble fully matured mRNAs with 5'cap1 structure and 3' polyA tail, therefore ready to be translated by the ribosome. mRNA transfection provides several advantages over plasmid DNA (pDNA) delivery. It does not require nuclear uptake for being expressed since translation of mRNA occurs directly into cytoplasm. Indeed, nuclear delivery (transport through nuclear membrane) is one the principal barriers for transfecting slow or non-dividing cells and consequently, mRNA transfection is particularly attractive for such purpose. This approach presents also the advantage of being non-integrative which is particularly appealing for stem cells, regenerative medicine or vaccine fields. Contrary to pDNA, mRNA cannot lead to genetic insertion causing mutations. Moreover, the protein expression from the mRNA is promoter-independent and faster than with DNA. For transfection we recommend RmesFect™ (#RM21000) and RmesFect™ Stem (#RS31000).

Luciferase detection

For transfections performed with MRNA20 or 21, the detection of Renilla Luciferase can be monitored using a bioluminescence assay kit.

Kit contents

R-Luc mRNAs-20: 20 µg of mRNA unmodified or modified.

R-Luc mRNAs-100: 100 µg of mRNA unmodified or modified.

R-Luc mRNAs-1000: 1 mg of mRNA unmodified or modified.

Storage

R-Luc mRNAs must be stored at -80°C. We recommend to aliquot the mRNA solution for a better storage.

OZ Biosciences SAS

163 avenue de Luminy
 Case 922, zone entreprise
 13288 Marseille cedex 09 - FRANCE
 Ph: +33 (0) 486 948 516

Fax: +33 (0) 463 740 015
contact@ozbiosciences.com
order@ozbiosciences.com

OZ Biosciences INC

7975 Dunbrook Road,
 Suite B
 San Diego CA 92126- USA
 Ph : + 1-858-246-7840
 Fax : + 1-855-631-0626

contactUSA@ozbiosciences.com
orderUSA@ozbiosciences.com

Related Products

Ref	Description
MRNA12	mRNA F-Luc (5moU)
MRNA12c	mRNA F-Luc (5moU)_Cy5
MRNA16	mRNA F-Luc unmod
MRNA16c	mRNA F-Luc unmod_Cy5
MRNA24	mRNA F-Luc (N1-mψ)
RM21000	RmesFect™ transfection reagent 1mL
RS31000	RmesFect™ Stem transfection reagent 1mL

Discover the complete list of mRNA at: www.ozbiosciences.com
Custom mRNAs are also available now!

Contact Us

Feel free to contact us for all complementary information and remember to visit our website to stay informed on the latest breakthrough technologies and updated on our complete product list. (www.ozbiosciences.com). For bulk, please contact us: order@ozbiosciences.com

Technical questions: tech@ozbiosciences.com
Order: order@ozbiosciences.com
General Information: contact@ozbiosciences.com

OZ Biosciences SAS

163 avenue de Luminy
Case 922, zone entreprise
13288 Marseille cedex 09 - FRANCE
Ph: +33 (0) 486 948 516
Fax: +33 (0) 463 740 015
contact@ozbiosciences.com
order@ozbiosciences.com

OZ Biosciences INC

7975 Dunbrook Road,
Suite B
San Diego CA 92126- USA
Ph : + 1-858-246-7840
Fax : + 1-855-631-0626
contactUSA@ozbiosciences.com
orderUSA@ozbiosciences.com