



CENTRIFUGES

# EBA 200

BLOOD / URINALYSIS / PEDIATRIC

## Press 'Start' in Case of Emergency

High-performance with built-in 8-place rotor, compact frame and all-metal containment for added safety. Ideal for blood, urine and pediatric tubes. STAT, PPP and coag. in minutes with the EBA 200 S model.

<b>MAX. CAPACITY</b>	8 x 15 mL
<b>MAX. RPM/RCF</b>	8,000/6,153



### QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



### FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



### FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



### FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



### ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.




### MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.

# HETTICH EBA 200 APPLICATION PACKAGES




## EBA 200 BLOOD TUBE PACKAGE - 8

	Tube Qty.	Volume (mL)	Size (mm)	200	200 S	200	200 S
				RPM		RCF	
	8	1.6 - 5	13 x 75	≤ 6,000	≤ 8,000	≤ 2,697	≤ 4,794
	8	4 - 7	13 x 100	≤ 6,000	≤ 8,000	≤ 3,461	≤ 6,153
	8	8.5 - 10	16 x 100	≤ 6,000	≤ 8,000	≤ 3,461	≤ 6,153
	4	8	16 x 125	≤ 6,000	≤ 8,000	≤ 3,461	≤ 6,153
Includes: (1) Centrifuge, (1) Built-in Angle Rotor, (8) Adapters 13 x 75 mm, (8) Adapters 13 x 100 mm							


PKG. NO. (EBA 200)	PKG. NO. (EBA 200 S)
200BLOOD8	200SBLOOD8

## EBA 200 PEDIATRIC TUBE PACKAGE - 8

	Tube Qty.	Volume (mL)	Size (mm)	200	200S	200	200 S
				RPM		RCF	
	8	0.5	10.7 x 36	≤ 6,000	≤ 8,000	≤ 2,214	≤ 3,935
	Includes: (1) Centrifuge, (1) Built-in Angle Rotor, (8) Adapters 0.5 mL						


PKG. NO. (EBA 200)	PKG. NO. (EBA 200 S)
200BLOOD8-PED	200SBLOOD8-PED

## EBA 200 URINALYSIS PACKAGE - 4

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
	4	12 KOVA	17 x 102	≤ 6,000	≤ 3,461
	4	15 Conical	17 x 120	≤ 6,000	≤ 3,461
Includes: (1) Centrifuge, (1) Built-in, Angle Rotor					

PKG. NO.
200URINALYSIS4

## EBA 200 CUSTOM PACKAGE

	<b>Supports over 19 application configurations, please inquire for more information</b>
	<b>Tube types include:</b> Standard, Blood, Pediatric, Urine, and Conical up to 15 mL



[hettweb.com/eba200](http://hettweb.com/eba200)

CONTACT YOUR LOCAL HETTICH REPRESENTATIVE TO LEARN MORE:



Discover more solutions on-line at [hettweb.com](http://hettweb.com) Toll Free: 1(866) 370-4388 Email: [info@hettweb.com](mailto:info@hettweb.com)  
© 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-EBA-200-E Rev A 01-2016.



LAB TECHNOLOGY