

# ecodyst

Accelerating the path to Discovery®

## EcoChyll X1 Condenser



Upgrades  
any existing  
rotovap

### Key Features

- Does not require dry ice or glycol
- Always available on-demand condenser
- Faster rates of evaporation
- Ready in 60 secs
- Only action required is to turn it ON

## Hydrogen Rotary Evaporator



All-in-one  
modern  
rotovap

- Built-in condenser offers required cooling temp
- Built-in vacuum pump controller
- 5L heating bath, room temp to 180°C
- Touchscreen display
- Motorized evaporating flask lift
- Speed 20–280 rpm, timer, interval operation
- Up to 3000mL evaporator flask size
- Chemical-resistant stainless steel spring PTFE vapor tub

**Our devices are more than twice  
as fast as traditional rotovaps.**

## Some of Our Valued Clients



Smart Self-cooling



Fast Evaporation



Small Footprint



Energy Efficient



Eco-friendly



United Kingdom:  
slsinfo@scientific-labs.com

International:  
slsexport@scientific-labs.com

Northern Ireland:  
NI.sales@scientific-labs.com

Republic of Ireland:  
ROI.sales@scientific-labs.com

