

ProSieve™ EX Stains


Revolutionary, Fast and Safe

These revolutionary stains provide faster staining times and less handling than other staining products. With special features that make each product unique, choosing the right stain for your research is easy.

- ProSieve™ EX Safe Stain; the ultimate fast solution that can provide a one step, safe stain in 10 minutes.

 2°C to 8°C

ProSieve™ EX Safe Stain with PAGEr™ EX Gels Offers Better Results in Half the Time



15 min gel run with PAGEr™ EX Gels and ProSieve™ EX Running Buffer at 275V, 10 minutes ProSieve™ EX Safe Stain, (total 25 minutes)

Ordering Information – ProSieve™ EX Stains

Cat. No. NA	Cat. No. EU	Product Name	Product Description	Size
223566	223566	ProSieve™ EX Safe Stain	One step, ten minute protein stain that is non-toxic	1 L
223567	223567	ProSieve™ EX Safe Stain	One step, ten minute protein stain that is non-toxic	25 mL