

Manual Tube Pickers & Manual Decappers

A manual device, with one-handed operation, the XTP-1 incorporates an eject button enabling tube pick and place operations whilst avoiding contact between the hand and tube.

The XTP-1 minimizes the risk of heat transfer to a frozen sample, as well as minimizing the risk of possible cross contamination.

Simple to use, even when wearing laboratory or cryogenic-handling gloves, the tube picker is available in either 96-rack or 48-rack format.



KEY FEATURES

- ✔ Single-handed operation, soft grip handle
- ✔ Easy to use whilst wearing gloves
- ✔ All components that potentially come into contact with samples are made from temperature-resistant and chemically-resistant plastics
- ✔ Tubes can be selected from any position in the rack
- ✔ Removes the need for manual contact with tubes reducing the risk of heat transfer to frozen samples and cross contamination
- ✔ Ideal for use with FluidX sample storage tubes and compatible with sample storage tubes from most other manufacturers

Specification

96-format XTP-1	length 150mm	diameter 17mm	weight 20g
48-format XTP-1	length 165mm	diameter 17mm	weight 23g

Manual Decappers: The ideal range of accessories for the manual removal of TPE caps and Screw caps.

Our TPE septum cap options comprise: A 1-way decapper for the individual removal of TPE caps and an 8-way decapper for the simultaneous removal of 8 TPE caps. For screw caps, we offer our FluidX screw cap Capper/Decapper, designed specifically for external thread tubes.

Ordering Information


10-5010	FluidX XTP-1 Manual Single Tube Picker, 96-format , 1 picker per case Suitable for FluidX and most other manufacturers 96-format tubes
10-5020	FluidX XTP-1 Manual Single Tube Picker, 48-format , 1 picker per case Suitable for FluidX and most other manufacturers 48-format Cryo Tubes
65-54000	FluidX Manual Decapper-1 , for removal of 1 TPE cap at a time, 1 decapper per case
65-54001	FluidX Manual Decapper-8 , for removal of 8 TPE caps at a time, 1 decapper per case
65-54004	FluidX Manual Screw Cap Capper/Decapper , 1 capper/decapper per case. Suitable for FluidX 96-format External Thread tubes



Automated Storage Systems



Cryopreservation & Cold Chain Solutions



Informatics & Technical Solutions



Sample Storage, Lab Services & Transport

Sample Consumables & Instruments