

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Specification

Product Name:

TEST

ANTI-FLAG* M2 Affinity Gel - purified immunoglobulin, buffered aqueous glycerol solution

Product Number: A2220

Appearance (Color)	White
Appearance (Form)	Suspension
Resin slurry in 10mM sodium phosphate, 150 mM sodium chloride,	pH 7.4, 50% (v/v) glycerol and 0.02% (w/v) sodium azide
Appearance	Pass
Upon settling: White resin with clear supernatant.	
Ligand Immobilized	> 2.7
grams antibody/liter gel	
Gel Content	1.7 - 2.1
Resin slurry will yield 1 ml PGV	
Binding Capacity / ml PGV	> 0.6
mg N-BAP eluted per ml of packed resin	_
Purity	> 95 %

Specification: PRD.1.ZQ5.10000020498

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.