

MSU-200 Series Compound Microscope

- Integrated X-Y mechanical stage
- Coaxial coarse and fine adjustment with rack stop
- Ergonomic stand
- Corded and cordless operation
- 5-year warranty



Cole-Parmer® MSU-200 Series Compound Microscope

Eyepiece(s)

- Wide field WF 10x/18 mm eyepieces

Head

- Binocular model with 30° inclined tubes, interpupillary distance is adjustable between 48 and 75 mm, and equipped with a ±5 diopter adjustment on the left tube
- 360° rotatable and secured with a screw

Nosepiece

- Reversed quadruple nosepiece on ball bearings

Objectives

- Semi-plan 4x/0.10, 10x/0.25, S40x/0.65 objectives, and S100x/1.25 oil immersion objective
- The 40x and S100x oil immersion objectives are spring loaded
- All optics are anti-fungus treated with anti-reflection coating for maximum light throughput

Stage

- 130 x 130 mm with mechanical 70 x 28 mm X-Y translation stage

Focusing

- Coaxial coarse and fine adjustments with 200 graduations, precision 2 µm per graduation, 0.4 mm per rotation
- Coarse adjustments equipped with friction control
- Total travel range approximately 23 mm
- Adjustable rack stop to prevent damage to samples and objectives

Condenser

- Height-adjustable Abbe condenser N.A. 1.25 with iris diaphragm and filter holder

Illumination

- Adjustable 1 W NeoLED illumination system for increased light output
- Rechargeable batteries and integrated power supply for cordless operation

What's Included

- Microscope, dust cover, spare fuse, 5 mL immersion oil for S100x objective, crosshairs eyepiece for centering, and power cord.

Description	SMP objectives	Objective type	Power	Dimensions	Weight	Item number
Binocular compound microscope	4/10/S40/S100x	Brightfield semi-plan	100 to 240 V, 50/60 Hz	14.2 x 6.9 x 8.9" (360 x 176 x 227 mm)	8.8 lb (4.0 kg)	78904-00

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

4932

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480.272279
All other countries: +1.847.549.7600

Cole-Parmer®