

KERN 317-01

Single weight, cylindrical, polished stainless steel

KERN[®]



Construction

Material:	polished stainless steel
Shape:	cylindrical

Measuring system

Nominal value:	1 g
OIML class:	E2
Tolerance:	0,030 mg

Accessories

Plastic box:	317-020-400
Wooden box:	317-010-100

Service

DAkkS Certificate:	962-331
Verification:	952-331

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	1 x 1 x 1 mm
Gross weight:	1 g
Net weight:	1 g
Shipping method:	Parcel service