

Rotors - fixed angle



BRK5401

max speed 12000 Rpm 13,500G
24 x Capillary tubes
Tube size maximum H75 x W1.6mm
Also 12 x 1.5-2.0ml tube
Tube size maximum H40 x W11mm
Cap size dependant.



BRK5424

Reducers available

max speed 12000 Rpm 13,500G
24 x 1.5-2.2ml tube, fixed angle rotor
Tube size maximum H40 x W11mm
Cap size dependant.
0.2ml - 0.4ml and 0.5ml



BRK5436

max speed 12000 Rpm 13,500G
36 x 0.5ml tube, fixed angle rotor
Tube size maximum H40 x W8mm
Cap size dependant

BRK5448

48 x 0.2ml tube, fixed angle rotor
Tube size maximum H40 x W6mm.



BRK5494

max speed 12000 Rpm 13,500G
4 x 8 place strip PCR rotor (0.2ml)



BRK5412

Reducers available

max speed 6000 Rpm 3,800G
12 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
5ml and 10ml



BRK5406

Reducers available

max speed 6000 Rpm 3,800G
6 x 50ml Tube, fixed angle rotor
Tube size maximum H115 x W30mm
15ml and 25ml

