


Serological pipette 2ml ClearLine® sterile in individual packaging**Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: green

- Sterile (Gamma)
- DNase and RNase free
- Non-pyrogenic

- Volume 2 ml
- Graduations: 0.02 ml
- Negative graduations 0.8 ml

Total length (A)	271 mm
Ø Tip (B)	2.7 mm
Ø Tip (C)	6.3 mm

**Conditioning**

- Individual paper/plastic packaging
- Box of 8 bags of 100 units

Units :

800 Tips