


**Serological pipette 5ml ClearLine® sterile in individual packaging****Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: blue
  
- Sterile (gamma)
- DNase and RNase free
- Non-pyrogenic
  
- Volume 5 ml
- Graduations: 0.1 ml
- Negative graduations 2.5 ml

Total length (A)	340 mm
Ø Tip (B)	3.2 mm
Ø Tip (C)	8.1 mm

**Conditioning**

- Individual paper/plastic packaging
- Box of 6 bags of 50 units

**Units :**

300 Tips