


Clearline tips 10 μ L Hybridrack ClearLine® sterile

- In clear polypropylene
- Store in a clean and dry place, at room temperature
- Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

| | |
|----------------------|---|
| Part number | 258361 |
| Packaging | HybridRack: cardboard rack with recyclable tip tray. (10 racks of 96 tips) |
| Volume | 10 μ l |
| Graduation | 2, 5, 10 |
| Length | 45,34 mm |
| Color | Natural |
| Filter | No |
| Low Retention | Yes |
| Sterility | Yes |
| Autoclavable | Yes |



Units :
960 Tips