


TECHNICAL DATASHEET**PART NUMBER****713246**

The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.

**ClearLine Tips**

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	713246
Packaging	ClearLine plus rack with hinged lid
Volume	10 µl
Graduation	1, 2, 5, 10
Length	44,8 mm
Color	Natural
Filter	No
Low Retention	No
Sterility	Yes
Autoclavable	Yes

Sales units :

960 tips