

TECHNICAL DATASHEET

PART NUMBER

713262




The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	713262
Packaging	ClearLine plus rack with hinged lid
Volume	100 µl
Graduation	No
Length	50,7 mm
Color	Natural
Filter	Yes
Low Retention	Yes
Sterility	Yes
Autoclavable	No



Sales units :

960 tips