

# COOLSAFE 55-9 SPECIFICATIONS



For more information please visit:  
[www.labogene.com](http://www.labogene.com)

CoolSafe 55-9	Unit	Part no 700100055
Digital read-out of temperature		Yes
Digital read-out of pressure		No
Ultimate temperature at 20°C room temperature	°C	-55
Power requirements	V/Hz	230/50-60
Power consumption	W	600
Insulation	mm	50
Condenser capacity	kg	6
Condenser capacity per 24 hours	kg	4
Total volume of Condenser	L	9
Weight net	kg	45
Cabinet dimensions (DxWxH)	mm	500x400x520
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507
Vacuum Valve		Solenoid
Drain tap		Yes
Keypad for control options		No
R 485/USB port for Labogene software		No

LaboGene ApS  
Bjarkesvej 5  
DK-3450 Allerød

Tel (+45) 3940 2566  
Fax (+45) 3995 2566  
Mail [info@labogene.com](mailto:info@labogene.com)  
Web [www.labogene.com](http://www.labogene.com)

