

# COOLSAFE TOUCH 100-9

## SPECIFICATIONS



For more information please visit:  
[www.labogene.com](http://www.labogene.com)

CoolSafe Touch 100-9	Unit	Part no: 700053
Digital read-out of temperature		Yes
Digital read-out of pressure		Yes
Ultimate temperature at 20°C room temperature	°C	-100
Power requirements	V/Hz	230/50-60
Power consumption	W	1300
Insulation	mm	90
Condenser capacity	kg	7
Condenser capacity per 24 hours	kg	4
Total volume of Condenser	L	9
Weight net	kg	75
Cabinet dimensions (DxWxH)	mm	620x495x745
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507/R170
Vacuum Valve		Electric
Drain tap		Yes
Touch screen for control options		Yes
USB port for LaboGene software/SD card for data		Yes

LaboGene ApS  
 Bjarkesvej 5  
 DK-3450 Allerød

Tel (+45) 3940 2566  
 Fax (+45) 3995 2566  
 Mail [info@labogene.com](mailto:info@labogene.com)  
 Web [www.labogene.com](http://www.labogene.com)

