

COOLSAFE TOUCH 55-15

SPECIFICATIONS



For more information please visit:
www.labogene.com

| CoolSafe Touch 55-15 | Unit | Part no: 700055 |
|---|------|--------------------------|
| Digital read-out of temperature | | Yes |
| Digital read-out of pressure | | Yes |
| Ultimate temperature at 20°C room temperature | °C | -55 |
| Power requirements | V/Hz | 230/50-60 |
| Power consumption | W | 800 |
| Insulation | mm | 90 |
| Condenser capacity | kg | 10 |
| Condenser capacity per 24 hours | kg | 6 |
| Total volume of Condenser | L | 15 |
| Weight net | kg | 70 |
| Cabinet dimensions (DxWxH) | mm | 620x495x745 |
| Cabinet material | | Polyester coated steel |
| Condenser material | | Stainless steel AISI 316 |
| Cooling media | | R 507 |
| Vacuum Valve | | Electric |
| Drain tap | | Yes |
| Touch screen for control options | | Yes |
| USB port for Labogene software/SD card for data | | Yes |

LaboGene ApS
 Bjarkesvej 5
 DK-3450 Allerød

Tel (+45) 3940 2566
 Fax (+45) 3995 2566
 Mail info@labogene.com
 Web www.labogene.com

