



designed for scientists



MMT 40.1 Multi-use milling tube

/// Data Sheet

Multi-use Milling Tube MMT 40.1

The Multi-use Milling Tube MMT40.1 with a maximum volume of 40 ml can be used and cleaned in a dishwasher multiple times. The package includes a rich set of spare parts so that wearing parts can be replaced if necessary. Therefore it is possible to carry out a number of experiments with the grinding chamber depending on the nature of the sample.

Delivery: 1 x MMT 40, 25 x sealing, 5 x beater, 5 x coupling

www.ika.com

Subject to technical changes





designed for scientists

Additional there is an Abrasion set A-MMT.100 available.



designed for scientists

Technical Data

| | |
|--------------------------------------|--------|
| Weight [kg] | 0.262 |
| Permissible ambient temperature [°C] | 5 - 40 |