

IKA

designed for scientists



MultiDrive MI 250.2 Star shaped cutter, stainless steel

/// Data Sheet

Star shaped cutter, stainless steel

Suitable for fibrous materials such as paper and plant material

Compatible with IKA MultiDrive MI 250 / 250 T



designed for scientists

Technical Data

Process type	batch
Operating principle	cutting
Feed hardness max. [Mohs]	5
Material beater/cutter	stainless steel 1.4034
FDA conform	yes
Multi use	yes
Beater/Knife diameter [mm]	66.5
Permissible ambient temperature [°C]	5 - 40