



STANDARD series

Underbench -20°C Freezer Model: BIO110FZSS

The **BioCOLD STANDARD** series of refrigeration has been built to a very high specification to meet the requirements for general laboratory storage.

The Underbench -20°C Freezer features 4 solid removable baskets enabling the safe storage of small vials. An optional inventory system for sample storage is available in various configurations.



Standard Features

- Sparkfree interior
- Reversible doors
- Replaceable door gasket
- Solid removable baskets
- CFC and HCFC free refrigeration system

Options

- Door mounted temperature display/alarm module
- Door lock

Technical Specification

Dimensions

External (H x W x D), (mm)
Weight (kg)

845 x 545 x 570
39

Temperature

Set Point (adjustable) (° C)
Operating range (° C)

-20
-18 to -23

Capacity

Freezer Gross, (L)

110

Electrical Supply

Voltage (V)
Frequency (Hz)
Fuse (A)

220-240
50
13

Drawers (solid baskets)

4

Defrost Type

Manual

