



# 24L

## Protected Access

Product Code:

### 5120XXXX

*Please note the last four digits denote the specific aperture type and lid colour.*

Our Sharpsafe containers are specifically designed to enable the safe disposal of clinical sharps and minimise the risk of associated injuries. Protected access offers extra peace of mind for patient safety, the additional protective shield prevents hands from entering the aperture. Ideally suited for high use areas and areas where children are present. The barrier also discourages non-sharps being disposed of in the containers, ensuring proper disposal and lowering associated costs.

Our 24L container is designed for high volume areas, such as A&E, wards and phlebotomy clinics. The wide horizontal aperture allows for the safe disposal of larger sharps equipment, surgical and laparoscopic instruments. Available with different colour lids for easily identifiable waste segregation. Products come with a range of compatible accessories for increased stability, security and safety.

Sharpsafe's new 5th Generation containers launched in 2022 and were the world's first full range of sustainable single use sharps containers. Engineered using ultra low carbon footprint recycled material showing our commitment to championing change.



SUPPORTING  
A CIRCULAR  
ECONOMY



REDUCING  
YOUR CARBON  
FOOTPRINT



REDUCING VIRGIN  
PLASTIC ENTERING  
THE SUPPLY CHAIN



MADE USING LOCALLY  
SOURCED LOW  
CARBON FOOTPRINT  
RECYCLED MATERIALS





Product Certification	
ISO 23907-1:2019	✓
UN3291	✓
BSI Kitemark	✓

Site Certification	
ISO 9001 - Quality Management	✓
ISO 14001 - Environmental Management	✓
ISO 45001 - Occupational Health & Safety	✓

Product Details	
Size:	24L
Usable Capacity:	18.85L
Usable Volume:	79%
Weight:	1.18kg
Fill Line Height:	236mm
Autoclavable:	Yes
Autoclavable Conditions:	134°C
Product Dimensions (mm):	335 x 450 x 230 (HxWxD)
Aperture Dimensions (mm):	365 x 95 (LxW)
GTIN:	5056271505572
Lot/Batch Number:	Not Applicable
Expiry Date:	Not Applicable

Accessory Compatibility			
Wall Bracket	✓	Pedestal Trolley	✓
Stand	n/a	Sharpsafe Trolley	✓
NPSD Tray	n/a	XL Trolley	n/a

Logistic Data	
Carton Quantity:	15
Carton Weight (when empty):	1.36kg
Carton Weight (when full):	17.7kg
Carton Dimensions (mm):	749 x 454 x 454 (HxWxD)
Carton Volume:	0.159m <sup>3</sup>
Storage Temperature Conditions:	Ambient
Storage Lower Temp Limit:	None
Storage Upper Temp Limit:	None
Sterile/Non Sterile:	Non Sterile

Features and Benefits			
In-Mould Label	✓	Latex Free	✓
Temporary Closure	✓	Low Carbon Footprint Material	✓
Permanent Closure	✓	Leak Resistance	✓
Four Click Assembly	✓	Disconnection Features	n/a
Fill Line	✓	Closure Shield	✓
Protected Access*	✓	Semi Translucent Lid	✓
Single Use	✓	Carry Handle	✓
Different Lid Variations**	✓	Puncture Resistant	✓

Pallet Information	
Euro Pallet Carton Quantity:	8
Euro Pallet Dimensions (mm):	1200 x 800
Standard Pallet Carton Quantity:	12
Standard Pallet Dimensions (mm):	1200 x 1000

Manufacturing Information	
Made In:	Wales, UK
By:	Vernacare

\* Additional product data sheet available  
 \*\* Different product codes may apply

FOR MORE INFORMATION:

✉ [sharpsafesales@vernagroup.com](mailto:sharpsafesales@vernagroup.com)

🌐 [vernacare.com](http://vernacare.com) 📞 +44 (0)1495 235800