



PRODUCT SPECIFICATION

TIE-ON FACE MASK

STERILE/NON-STERILE

FaceGuard **1**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

A 3 layer tie on facemask with flexible nose piece to ensure an improved fit and wearer comfort

PRODUCT FEATURES

Outer layers manufactured from Ethylene-Propylene Side by Side (ES)
Middle layer manufactured from melt blown Polypropylene
Flexible nose-piece for improved fit and wearer comfort
Polyester Tie-on loop provided to give improved adjustment
Size: 21 x 9 cm

MANUFACTURING

Manufactured and packed in a cleanroom environment

PACKAGING

Sterile variant individually packed and double-bagged

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy) sterile variant only

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems

SHELF LIFE

2 years from date of manufacture

STORAGE CONDITIONS

Store in cool, well ventilated area

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
	MSTW01020	Individually Packed Tie-on Face Mask Sterile	21 x 9cm	20	15 x 20	310 x 270 x 373mm	2.4kg
MNTW01100		Multi Packed Tie-on Face Mask Non-Sterile	21 x 9cm	100	10 x 100	520 x 270 x 350mm	6.7kg
MNTW01020		Individually Packed Tie-on Face Mask Non-Sterile	21 x 9cm	20	50 X 20	520 x 560 x 350mm	7.8kg