

|   |                                |                      |           |
|---|--------------------------------|----------------------|-----------|
|  | <h1>Product Specification</h1> | Specification Number | 8045A-500 |
|   |                                | Revision             | 07        |

Part No. 8045A-500

Description SteriWare PharmaScoop® – 500ml

Use: Single use, disposable scoop

Material of Construction: High Impact Polystyrene (HIPS) – rigid, white

- FDA acceptable CFR 177.1640
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant



Method of Construction: Single piece, injection moulded

Overall Length: 249mm

Handle Length: 118mm

Maximum Width: 97mm

Nominal Weight: 86g

Assembly & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: N/A

Number of Samplers per Box: 25 (Supplied as 4 boxes of 25 pcs)

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare and PharmaScoop are registered Trademarks belonging to Sampling Systems Ltd.

|                            |  |   |
|----------------------------|--|---|
| Specification Approved By: | Signature:<br> | Name: Jonathan Fenner<br>Position: Director<br>Date: 3-Oct-17 |
|----------------------------|--|---|