

1. Part Number: 8088A-03S

2. Photograph:



3. Product Information:

Description:	SteriWare® Sample Spoon 3ml - Sterile
Use:	Single use, disposable spoon
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1640 • Conforms to EU regulation 10/2011 • Conforms to EC regulation 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Bowl Volume:	3ml
Length:	173mm
Bowl Length:	36mm
Bowl Width:	25mm
Nominal Weight:	4g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 19-Jul-19