



1. Part Number: 1220J-2000

2. Photo:



3. Product Information:

Description	Jumbo Pocket Sampler 2000mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Internal Polish:	None
Overall Length:	2050mm
Diameter of Outer Tube:	40mm
Pocket Length:	125mm
Pocket Width:	20mm
Nominal Pocket Volume:	90ml
Number of Pockets:	4
Slot Positions:	



Nominal Weight:	7.3 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Mar-20