

**1. Part Number:** 8045R-60

**2. Photograph:**



**3. Product Information:**

|   |  |
|---|--|
| Description:                            | SteriWare PharmaScoop® RED 60ml  |
| Use:                                    | Single use, disposable scoop   |
| Material of Construction:               | High Impact Polystyrene (HIPS) – rigid, red <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1640</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC regulation 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:                 | Single piece, injection moulded  |
| Nominal Volume:                         | 60ml   |
| Overall Length:                         | 178mm  |
| Handle Length:                          | 106mm  |
| Maximum Width:                          | 49mm   |
| Nominal Weight:                         | 14g  |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom  |
| Individually Bagged?                    | Yes (heat sealed PE/PET bag)   |
| Method of Sterilisation:                | N/A  |
| Number of Scoops per Box:               | 100  |
| Recommended Storage Conditions:         | Dry and ambient temperature  |
| Shelf Life:                             | 5 years from the date of manufacture   |
| ATEX:                                   | This device has not been ATEX rated  |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Dec-19