

1. Part Number: A716-100

2. Photograph:



3. Product Information:

Description:	Powder Funnel 100mm
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	340ml
Top Diameter:	100mm
Overall Height:	107mm
Spout Length:	25mm
Internal Spout Diameter:	35mm (Spout ends with a 1.5" tri-clamp ferrule BS4825-3)
Nominal Weight:	460g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Aug-20