

1. Part Number: A716-400-2

2. Photograph:



3. Product Information:

Description:	Powder Funnel 400mm with 2" ferrule
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	23.3 Litres
Top Diameter:	400mm
Top Diameter including flange:	418mm
Overall Height:	396mm
Spout Length:	25mm
Spout Diameter:	Spout ends with a 2" ferrule (BS4825-3)
Nominal Weight:	4.2kg
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Jul-20