

|   |                                |                      |           |
|---|--------------------------------|----------------------|-----------|
|  | <h1>Product Specification</h1> | Specification Number | 8095B-250 |
|   |                                | Revision             | 01        |

Part No. 8095B-250

Description Disposable Powder Spatula Blue – 250mm

Material of Construction: Polypropylene (PP) – rigid, blue

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004



Method of Construction: Single piece, injection moulded

Nominal Length: 330mm

Blade Length: 246mm

Blade Width: 18mm

Nominal Weight: 22g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: N/A

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade marks belonging to Sampling Systems Ltd.

|                            |  |  |
|----------------------------|--|--|
| Specification Approved By: | Signature:<br>   | Name: Jonathan Fenner<br>Position: Director<br>Date: 21-Apr-17 |
| Page 1 of 1                | <b>Sampling Systems Ltd.</b><br>4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK.<br>Tel: +44(0)1675 466992 Fax: +44(0)1675 466994<br>Email: info@sampling.com Web: www.sampling.com |  |