

**1. Part Number:** 8064H-1000

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® Powder Probe 1000mm
Use:	Single use powder sampler
Material of Construction:	High Density Polyethylene (HDPE) – virgin <ul style="list-style-type: none"><li>• Conforms to FDA 21 CFR 177.1520</li><li>• Conforms to EU Regulations 10/2011</li><li>• Conforms to EC Regulations 1935/2004</li><li>• BSE/TSE free</li></ul>
Method of Construction:	Ultrasonic welding
Nominal Length:	945mm
Blade Length:	845mm
Diameter of Blade:	22mm
Width of opening in blade:	10mm
Nominal Weight of Sampler:	123g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

**4. Change Control:** Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:** Made and packed in the UK.

## 6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

## 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21