

1. Part Number: 8284R-250

2. Photograph:



3. Product Information:

| | |
|---|--|
| Description: | SteriWare® Cupped Spatula Red – 250mm |
| Use: | Single use, powder spatula |
| Material of Construction: | Polystyrene – virgin, red <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1640 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Blade Shape: | U |
| Blade Style: | Closed Tip |
| Overall Length: | 335mm |
| Blade Length: | 250mm |
| Blade Width: | 18mm |
| Nominal Weight: | 24g |
| Nominal Volume (flat-filled): | 40ml |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | N/A |
| Number of Samplers per Box: | 100 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This item has not been ATEX tested |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20