

**1. Part Number:** 8286A-1000S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Scoop 1000ml - Sterile
Use:	Single use, disposable scoop
Nominal Flat-filled Volume:	1000ml
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1640</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Overall Length	360mm
Handle Length:	134mm
Maximum Width:	127mm
Nominal Weight:	119g
Moulding & Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Scoops per Box:	30
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable scoop. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Jul-19