

HyClone classical media

CELL CULTURE

Dulbecco's modified eagle's medium - low glucose Liquid media

Component description	Code number SH30021, SH30307	
	mg/L	mmol/L
Inorganic salts		
Calcium chloride	200	1.8021
Ferric nitrate-9H ₂ O	0.1	0.0002
Potassium chloride	400	5.3655
Magnesium sulfate	97.67	0.8112
Sodium chloride	6400	109.514
Sodium phosphate monobasic H ₂ O	125	0.9059
Amino acids		
L-Arginine-HCl	84	0.3987
L-Cystine-2HCl	62.57	0.1998
L-Glutamine	584	3.9959
Glycine	30	0.3996
L-Histidine-HCl-H ₂ O	42	0.2004
L-Isoleucine	104.8	0.7989
L-Leucine	104.8	0.7990
L-Lysine-HCl	146.2	0.8004
L-Methionine	30	0.2011
L-Phenylalanine	66	0.3995
L-Serine	42	0.3997
L-Threonine	95.2	0.7992
L-Tryptophan	16	0.0783
L-Tyrosine-2Na-2H ₂ O	103.79	0.3974
L-Valine	93.6	0.7990

Component description	Code number SH30021, SH30307	
	mg/L	mmol/L
Vitamins		
Calcium D-pantothenate	4	0.0084
D-Pantothenic acid Na salt	0	0
Choline chloride	4	0.0286
Folic acid	4	0.0091
Myo-inositol	7	0.0389
Niacinamide	4	0.0328
Pyridoxine-HCl	4	0.0195
Riboflavin	0.4	0.0011
Thiamine-HCl	4	0.0119
Other		
D-Glucose	1000	5.5506
HEPES	0	0
Phenol red-Na	15.9	0.0422
Sodium pyruvate	110	0.9996
Sodium bicarbonate	3700	44.0424

[cytiva.com/hyclone](https://www.cytiva.com/hyclone)

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. HyClone is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY13873-04Aug20-DF

