



designed for scientists



AS 501.5 Dish attachment

/// Data Sheet

For smoothly shaking dishes, but also for smooth mixing in vessels with a large, flat bottom (widenecked Erlenmeyer flasks and beakers). A plastic foil with mild adhesive prevents the vessel from slipping.

Dimensions (W x D x H): 450 x 450 x 45 mm

Set-up plate: 420 x 420 mm

Weight: 1.120 g



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	505 x 480 x 225
Weight [kg]	0.987
Permissible ambient temperature [°C]	5 - 40