



AS 4000.3 Dish attachment

/// Data Sheet

For smooth shaking operations in the low viscosity range, e.g. for cell cultures, nutrient media in Petri dishes, culture bottles and vessels with a low center of gravity. With integrated slip-resistant foil (PP).

Dimensions (W x D x H): 470 x 444 x 25 mm

Set-up plate: 430 x 430 mm

Weight: 800 g

Clamping screws for fixing the attachment are included in the scope of delivery of the KS 4000 incubator shaker.



designed for scientists

Technical Data

Depth [mm]	150
Height [mm]	60
Weight [kg]	1.39
Permissible ambient temperature [°C]	5 - 40