



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: ACETIC ANHYDRIDE A.R.

CAS No: 108-24-7

<i>Supplied as:</i>	Product code	Description
	CHE1024	ACETIC ANHYDRIDE A.R. 2.5L.
	CHE1026	ACETIC ANHYDRIDE A.R. 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98.5			
Wt. per ml.	g	1.079	1.082		
Chloride	%		0.0002		
Sulphate	%		0.0005		
Non-volatiles	%		0.0050		

Revision date: 27/10/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0013/MS/20211027