



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: ACETONE pure

CAS No: 67-64-1

<i>Supplied as:</i>	Product code	Description
	CHE1040	ACETONE pure 2.5L.
	CHE1042	ACETONE pure 25L.
	CHE1044	ACETONE pure 1L.
	CHE5084	ACETONE pure 5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.0			
Wt. per ml.	g			0.790	

Revision date: 13/07/1990

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0014/MS/19900713