

Material description: AMMONIUM CHLORIDE A.R.

CAS No: 12125-02-9

Supplied as: Product code Description
CHE1150 AMMONIUM CHLORIDE A.R. 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
pH; 5% solution	pH	4.5	5.5		
Calcium	%		0.0005		
Iron	%		0.0002		
Phosphate	%		0.0002		
Sulphate	%		0.002		
Insoluble matter	%		0.005		
Heavy Metals	ppm		5		
Magnesium	ppm		5		

Revision date: 24/06/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3381/MS/20210624