



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: CHOLESTEROL

CAS No: 57-88-5

Supplied as: Product code Description  
CHE1580 CHOLESTEROL 25g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Cholesterol content	%	90.0			
Water content <i>by Karl Fischer</i>	%		0.5		
Appearance					White / Off-white powder

Revision date: 07/07/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0178/MS/20210707