



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: 1-IODOBUTANE pure

CAS No: 542-69-8

<i>Supplied as:</i>	Product code	Description
	CHE2226	1-IODOBUTANE pure 25ml.
	CHE2228	1-IODOBUTANE pure 50ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Boiling range	°C	129	132		
Wt. per ml.	g	1.60	1.64		

Revision date: 23/07/1990

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0378/MS/19900723