

Material description: PROPAN-2-OL pure  
CAS No: 67-63-0

<i>Supplied as:</i>	Product code	Description
	CHE3116	PROPAN-2-OL pure 2.5L.
	CHE3118	PROPAN-2-OL pure 500ml.
	CHE3120	PROPAN-2-OL pure 25L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Boiling range	°C	81	83		
Wt. per ml.	g	0.784	0.787		

Revision date: 04/02/1994

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0557/MS/19940204