

Material description: SAND PURIFIED (dried white silica sand)  
CAS No: 7631-86-9

<i>Supplied as:</i>	Product code	Description
	CHE3184	SAND PURIFIED (dried white silica sand) 1Kg.
	CHE3186	SAND PURIFIED (dried white silica sand) 5Kg.
	CHE3188	SAND PURIFIED (dried white silica sand) 25Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Silica	%			99.7	
Aluminium Oxide	%			0.07	
Ferric Oxide	%			0.03	
Sodium Oxide	%		0.05		
Potassium Oxide	%		0.01		
Calcium Oxide	%		0.02		
Magnesium Oxide	%		0.05		
Acid washed					

Revision date: 26/01/2012 Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1788/MS/20120126