

Material description: SILVER NITRATE A.R.

CAS No: 7761-88-8

<i>Supplied as:</i>	Product code	Description
	CHE3224	SILVER NITRATE A.R. 25g.
	CHE3226	SILVER NITRATE A.R. 100g.
	CHE3228	SILVER NITRATE A.R. 250g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.9			
Chloride	%		0.0005		
Nitrite	%		0.0005		
Sulphur	ppm		8.0000		
Copper	%		0.0005		
Iron	%		0.0005		
Lead	%		0.0005		

Revision date: 27/07/1993

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0578/MS/19930727