

Material description: TRIS-(HYDROXYMETHYL)-METHYLAMINE A.R.

CAS No: 77-86-1

<i>Supplied as:</i>	Product code	Description
	CHE3834	TRIS-(HYDROXYMETHYL)-METHYLAMINE A.R. 250g.
	CHE3836	TRIS-(HYDROXYMETHYL)-METHYLAMINE A.R. 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Melting range	°C	168	172		
Appearance					White crystals
pH; 5% solution	pH	10.0	11.5		
Chloride	%		0.001		
Sulphated ash	%		0.10		
Water	%		0.5		

Revision date: 12/10/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2407/MS/20201012