

Material description: tri-SODIUM CITRATE 2H₂O ACS
CAS No: 6132-04-3

Supplied as: Product code Description
CHE3838 tri-SODIUM CITRATE 2H₂O ACS 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Assay	%	99.00			
Heavy Metals	ppm		5		
Calcium	%		0.0030		
Chloride	%		0.003		
Iron	ppm		5		
Ammonia	%		0.003		
Sulphate	%		0.005		
Insoluble matter	%		0.005		
pH 5% aqueous solution	pH	7.00	9.00		

Revision date: 07/01/2022

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 4918/MS/20220107