

Material description: PROPANE-1,2-DIOL pure
CAS No: 57-55-6

<i>Supplied as:</i>	Product code	Description
	CHE3126	PROPANE-1,2-DIOL pure 500ml.
	CHE3128	PROPANE-1,2-DIOL pure 2.5L.
	CHE4040	PROPANE-1,2-DIOL pure 1L.
	CHE5072	PROPANE-1,2-DIOL pure 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Boiling range	°C	185	189		
Wt. per ml.	g	1.035	1.037		
Acidity; acetic	%		0.002		
Water	% w/w		0.2		

Revision date: 18/02/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1785/MS/20210218